



Septicare®

# HOOF GLUE SEPTICARE II BOND EXPRESS



**SINGLE-USE SYRINGE,  
PERFECT FOR 1 HOOF!**

## SEPTICARE II BOND EXPRESS

Sku. 106652005

Unchanged formula, same efficiency, more practicality!  
Applied in **3 minutes**, Stick in **30 seconds**,  
the cow is released in 3 minutes.



- **Significant time and cost savings:**  
the glue is applied in only 3 seconds,  
bonds in 30 seconds thanks to a very  
quick thermic reaction between the  
2 components. The cow is thus  
released in only 3 minutes.

- **A glue with unequalled efficiency:**  
the glue dries even under extreme  
conditions. It can be applied and  
stored from 0°C to 40°C and even  
in wet environments, but it must  
not have a temperature lower than  
0°C for application.

- **Very easy to use:**  
no need to be preheated,  
the cannula makes the 2 component  
optimal mixture for a maximal  
bonding and the glue doesn't  
flow even by applying vertically.

- **Comfort for the breeder and the cow:**  
the breeder saves 15 seconds  
per treated cow, the cow is less  
stressed and the treatment has  
no impact on the milk production.

**AVAILABLE IN 250ML  
SEPTICARE II BOND  
EXPRESS XL**

Sku. 106652001



## BOX OF 10 KITS

Each kit contains one 25 ml syringe + 1 cannula  
+ 1 pair of gloves, individually wrapped



3 700393 919718

**www.septicare.eu**

Order on our Webshop

→ **www.ukal-elevage.com**